Instructions for Use

Gebrauchsinformation

Mode d'emploi

Istruzioni d'uso

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- Light-curing resin-based dental restorative material
- Lichthärtender zahnärztlicher Füllungskunststoff
- Matériau de restauration en résine composite, photopolymérisable
- Materiale da restauro fotopolimerizzabile su base resinosa
- Material de restauración dental basado en resina fotopolimerizable
- Compósito fotopolimerizável para restauração dental

Bruksanvisning

Brugsanvisning

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Bruksanvisning

Productinformatie

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For dental use only.

Caution: U.S. Federal Law restricts this device to sale by or on the order of a licensed dentist.

Made in Liechtenstein

Ivoclar Vivadent AG, FL-9494 Schaan/Liechtenstein

Complies with / entspricht:

ISO 4049, EN ISO 4049

CE 0123

573263/0505/WE3/G

ivoclar
vivadent®
clinical

Instructions for Use

Description

Thank you for purchasing Artemis restorative material which enables you to meet the highest expectations in restorative therapy.

The aesthetic challenge of modern restoration techniques is to produce true-to-nature teeth. Artemis offers you a complex shade system and high-quality materials. Thus, the restorations imitate the optical properties of natural teeth more closely than ever.

Artemis is a radiopaque fine-particle hybrid composite for highly aesthetic restorative therapy. The material cures with light at a wavelength range of 400–500 nm (blue light).

Shades

Various dentin and enamel shades with different degrees of transparency are indispensable for highly aesthetic restorations. In combination with the available Effect shades, it is possible to imitate the natural anatomy in a nearly perfect way.

The following shades and degrees of transparency are available:

Dentin shades (Dentin)

A2, A3, A3.5, A4

B3

C4

D2, D4

IVA5, IVA6

Enamel shades (Enamel)

A1, A2, A3, A3.5, A4

B1, B2, B3, B4

C2

D2, D3

Bleach XL, L, M

Transparency

7–8 %

13–15 %

10–20 %

Effect shades (Effect)

White	6 %
Blue, Amber	21–26 %
Clear	30 %
Super Clear	50 %

Dentin, Enamel and Effect shades can easily be distinguished due to their optical appearance. Each group of shades is identified by a specific blue colour of the delivery form. Single shades can therefore be allocated to the three groups, which makes the system clear and user-friendly.

Composition

The monomer matrix is composed of Bis-GMA, urethane dimethacrylate and triethylene glycol dimethacrylate (22–24 wt%). The inorganic fillers contain barium glass, ytterbium trifluoride, Ba-Al-fluorsilicate glass, highly dispersed silicon dioxide and spheroid mixed oxide (75–77 wt%).

Additional contents: catalysts, stabilizers and pigments (<0.5 wt%). The total content of inorganic fillers is 75–77 wt% or 55–58 vol%.

The particle size is 0.04–3.0 µm. The mean particle size is 0.6 µm.

Indication

- Anterior restorations (Classes III, IV)
- Class V restorations
- Restorations in the posterior region (Classes I and II)
- Direct Veneering
- Diastema closure
- Correction of tooth position and tooth shape (eg closure of interdental triangles, lengthening of incisal edges)

Contraindication

Placement of Artemis is contraindicated:

- if a dry working field cannot be established, or if the stipulated technique cannot be applied
- if a patient is known to be allergic to any of the ingredients in Artemis

Side effects

In individual cases, components of Artemis may lead to sensitization in predisposed persons. Artemis should not be used in such cases.

To avoid possible irritation of the pulp, areas close to the pulp should be protected with a suitable pulp/dentin protector (see instructions under "Application" Point 4).

Interactions

Substances containing eugenol/oil of cloves may inhibit the polymerization of composite materials. Consequently, application of such materials together with Artemis must be avoided.

Discolouration may occur in combination with cationic mouthwashes, plaque disclosing agents, and chlorhexidine.

Application

1. Shade determination

Clean the teeth prior to shade determination with a polishing paste. The shade is selected with the tooth still moist.

The Artemis shade guide is useful for this purpose. The necessary shades for the dentin as well as the enamel can be selected in order to design and imitate the natural anatomy and optic appearance of the tooth.

On the Artemis shade guide, Dentin shades can be distinguished from the Enamel and Effect shades by the shape of the shade tab.

The colour code of the three groups in the delivery form applies to the shade samples too.

2. Isolation

Appropriate isolation, best with a rubber dam, is required.

3. Cavity preparation

Cavity preparation is carried out according to the requirements of the adhesive technique, i.e. protecting the tooth structure. Do not prepare sharp, internal edges or any additional undercuts in caries-free areas. The dimensions of the cavity are generally determined by the extent of the caries or the size of the old restoration. Bevel enamel edges in the anterior region to improve the aesthetic integration and marginal adaptation of the restoration.

In the posterior region, only the sharp enamel edges should be rounded (finishing diamonds, 25–40 µm). Caries-free cervical defects are not prepared, only cleaned with pumice or other suitable cleaning pastes with the help of rubber cups or rotary brushes. Subsequently, remove all residue in the cavity with water spray.

4. Pulp protection / Base

Do not apply a base material when applying an enamel/dentin bonding agent. Only cover very deep areas close to the pulp with a calcium hydroxide material and subsequently use a pressure-resistant cement (e.g. glass ionomer cement, such as Vivaglass Liner or zinc phosphate cement, such as PhosphaCEM PL).

Do not cover other cavity walls, since they can be used to support the bond with an enamel/dentin adhesive.

5. Matrix / Interdental wedge

Use a transparent matrix for cavities affecting the proximal area (e.g. Contour Strip in the anterior region) and wedge it. The use of steel matrices in the posterior region is also possible.

6. Conditioning / Application of the bonding agent

Conditioning and application of the bonding agent according to the Instructions for Use of the product in use. Given the coordination of Ivoclar Vivadent materials, we recommend using Syntac (with phosphoric acid etching) or Excite (with phosphoric acid etching) or the AdheSE self-etching adhesive.

7. Application of Artemis

As Artemis is available in a wide range of dentin and enamel shades and various degrees of transparency as well as in shades for creating anatomic details and special effects, any layering technique can be used to fabricate highly aesthetic restorations.

Apply Artemis in layers of max. 2 mm or 1.5 mm (dentin shades) and adapt with suitable instruments (e.g. P1 instruments). Sufficient exposure to the curing light prevents incomplete polymerization. Therefore, light cure

each increment with a polymerization unit with a light intensity of at least 500 mW/cm² (e.g. bluephase C5, C8) for 20 seconds. With a light intensity of at least 1100 mW/cm² (e.g. bluephase), the individual increments can be cured within 10 seconds.

Hold the light emission window as closely as possible to the surface of the restorative material.

When using a metal matrix, additionally polymerize from the buccal or lingual/palatal after removing the matrix.

In most cases, a flowable composite is used as first layer in order to create an even cavity bottom and to increase the primary bonding properties of the subsequently applied filling material. The application of a flowable composite (e.g. Tetric Flow) as a thin first layer is optional. This layer has to be polymerized separately. (Please refer to the respective Instructions for Use.)

8. Finishing / Checking the occlusion / Polishing

Remove excess material with suitable finishers (e.g. Astropol F) or fine diamonds after polymerization. Remove proximal excess with diamond or tungsten carbide finishers or finishing strips. Check occlusion and articulation and apply appropriate grinding corrections to prevent premature contacts or undesired articulation paths on the surface of the restorations. Use silicone (e.g. Astropol P, Astropol HP or Astrobrush), polishing discs, and polishing strips to polish the restorations to a durable high gloss.

Additional information

1. If necessary, additional Artemis can be directly applied to polymerized material. If the Artemis restoration has already been polished, it must first be roughened and wetted with Heliobond before fresh Artemis is applied.
2. Artemis should have room temperature when applied. Cool temperatures render the material difficult to extrude.
3. If Artemis from the Cavifill is applied directly in the mouth of the patient, we recommend using the Cavifill only once because of hygienic reasons (to avoid cross infections between patients).

Warning

Unpolymerized Artemis should not come in contact with the skin, mucous membrane, or eyes. Unpolymerized Artemis can have a slight irritating effect and may lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Storage

- Do not use Artemis after the indicated date of expiration.
- Store packages at 2–28 °C / 36–82°F.
- Close syringes / Cavifils immediately after use. Exposure to light causes premature polymerization.
- Shelf life: see information on syringes / packages.

Keep material out of children's reach.

For use in dentistry only.

Date information prepared

04/2005

Manufacturer

Ivoclar Vivadent AG

FL-9494 Schaan/Liechtenstein

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damage resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent UK Limited
Ivoclar Vivadent, Inc.
Dentigenix Inc.